REF RCUM-00171 / RCUM-00271 / RCUM-02071

ENGLISH

Quick Reference Instructions (QRI)

CareSuperb™ COVID-19/Flu A&B Antigen Combo Home Test

For *in vitro* diagnostic use
Over-the-Counter (OTC)
Please read this QRI carefully before use

Warnings and Precautions

- Do not use this test when you have symptoms greater for 4 days or more, or if you have no symptoms.
- Read all the instructions carefully before performing the test. Failure to follow the instructions may result in inaccurate test results.
- This test is read visually. Individuals with impaired vision or color-impaired vision should ensure help in interpretation of their test results.
- Testing should be performed in an area with good lighting.
- Taking biotin supplements may negatively affect the performance of this test. If you are taking biotin supplements and test negative, you should follow up with your physician.
- Keep the test cassette on a flat surface during the testing.
- Use the test components immediately once opened. Test components are single use. Do not re-use.
- Do not use if any of the test kit contents or packaging is damaged or open.
- When collecting a nasal swab sample, use ONLY the nasal swab provided in the kit.
- Do not use this test for individuals who recently received nasally administered influenza A or influenza B vaccine, as they may have false positive test results after vaccination.

- Keep foreign substances and household cleaning products away from the test during the testing process. Contact with foreign substances and household cleaning products may result in incorrect test results.
- Do not conduct this test if prone to nose bleeds or have a nose injury.
- Remove any piercings from nose before starting the test.
- In case of spillage, ensure it is cleaned thoroughly with a suitable disinfectant.
- Keep testing kit and kit components away from children and pets before and after use.
- Avoid contact with your skin or eyes. Do not ingest any kit components. The reagent solution contains harmful chemicals (see table below). If the solution contacts your skin or eyes, flush with large amounts of water. If irritation persists, seek medical advice: <https://www.poisonhelp.org> or 1-800-222-1222.

Chemical Name / Concentration	GHS Code / Harm
Triton x-100 / 1.0% Trizma base / 1.2% Sodium chloride / 2.9%	H315 / Skin irritation H320 / Eye irritation


- For the most up-to-date information on COVID-19, please visit: www.cdc.gov/COVID19.

PART A Prepare for the Test

1. Wash your hands with soap and water for at least 20 seconds and dry them thoroughly or use hand sanitizer.

2. Open your test kit and locate the kit components from the box.


3. Open the pouch and remove the test cassette from the foil pouch. Place the test device on a flat surface.


PART B Collect Nasal Sample


4. Remove a swab from the wrapper.

5. Gently insert the swab no more than 3/4 inch into the LEFT nostril. Then, slowly rotate the swab at least 5 times in a circular path for a total of 15 seconds.

6. Gently remove the swab from the LEFT nostril and place it directly into RIGHT nostril, repeating the process of rotating at least 5 times in a circular path for a total of 15 seconds. Remove the swab from the RIGHT nostril.


- Be careful not to touch the swab tip.

- Note: If you are swabbing others, please wear a face mask. With children, the maximum depth of insertion into the nostril may be less than 1/2 to 3/4 of an inch, and you may require another adult to hold the child's head while swabbing.


PART C Perform the Test

7. Insert the collected swab into the sample port and gently push the swab through the sealed silver foil sheet. Continue to gently push until the swab tip is FULLY NOT VISIBLE, ensuring it is completely below the sample port.


- Insert the swab into the sample port before using the dropper vial and keep the swab in the sample port until further instructions.

8. Remove the dropper vial cap by twisting it off while it's facing upward.

- Do not use mouth or teeth to open the bottle.

9. Add the ENTIRE VOLUME of the liquid into the sample port.

- Do not put the liquid in your eyes. Invalid or inaccurate results may occur if the entire volume is not added to the sample port. Make sure not to spill the liquid outside the sample port.

REF RCUM-00171 / RCUM-00271 / RCUM-02071

CareSuperb™ COVID-19/Flu A&B Antigen Combo Home Test

A rapid test for the detection of SARS-CoV-2 nucleocapsid antigens and Influenza A and B nucleoprotein antigens in anterior nasal specimens.

For *in vitro* diagnostic use
Over-the-Counter (OTC)

Intended Use

The CareSuperb™ COVID-19/Flu A&B Antigen Combo Home Test is a lateral flow immunochromatographic assay intended for the qualitative detection and differentiation of Influenza A and Influenza B nucleoprotein antigens and SARS-CoV-2 nucleocapsid antigen directly in anterior nasal swab samples from individuals with signs and symptoms of respiratory tract infection. Symptoms of respiratory infections due to SARS-CoV-2 and Influenza can be similar. This test is for non-prescription home use by individuals aged 14 years or older testing themselves, or adults testing individuals aged 2 years or older.

All negative results are presumptive and should be confirmed with an FDA-cleared molecular assay when determined to be appropriate by a healthcare provider. Negative results do not rule out infection with Influenza, SARS-CoV-2, or other pathogens.

Individuals who test negative and experience continued or worsening respiratory symptoms, such as fever, cough, and/or shortness of breath, should seek follow up care from their healthcare provider.

Positive results do not rule out co-infection with other respiratory pathogens and therefore do not substitute for a visit to a healthcare provider for appropriate follow-up.

Kit Components

Component Name	Quantity (in a kit)				Description
	Product Reference No. (REF)	RCUM-00171	RCUM-00271	RCUM-02071	
Test Device	1	2	20		Individually foil pouch test cassette
Dropper Vial	1	2	20		The vial containing extraction reagent
Nasal Swab	1	2	20		Sterile swab for anterior nasal specimen collection
Quick Reference Instructions (QRI)		1			Paper-printed instruction

*Materials required but not provided: Timer or watch

- Carefully read these QRI instructions before performing the test. Failure to follow the instructions may result in inaccurate results.
- A full IFU with performance data is available on the website at: www.accessbio.net/products/rapid-diagnostic-tests/caresuperb-covid-19-and-flu-a-and-b-antigen-combo-home-test.


Storage and Handling

- Store kit at 35.6 - 86°F (2-30°C).
- The expiration date is on the package. Do not use kit past its expiration date.
- The test cassette must remain in the sealed foil pouch until use. Once the pouch has been opened, it should be used immediately.
- Operate your test in an ambient temperature of 59°F to 86°F (15-30°C).

Limitations


- The performance of this test was established based on the evaluation of a limited number of clinical specimens collected between November 2023 and March 2025. The clinical performance has not been established for all circulating variants but is anticipated to be reflective of the prevalence variants in circulation at the time and location of the clinical evaluation. Performance at the time of testing may vary depending on the variants circulating, including newly emerging strains of SARS-CoV-2 and their prevalence, which change over time.
- There is a higher chance of false negative results with antigen tests than with laboratory-based molecular tests due to the sensitivity of the test technology. This means that there is a higher chance this test will give a false negative result in an individual with COVID-19, and Influenza as compared to a molecular test, especially in samples with low viral load.
- No pharmacokinetic studies are available to demonstrate that biotin may be secreted in respiratory samples after ingestion of high-dose biotin supplements however, studies have shown that false negative results occur when this test is exposed to concentrations of biotin exceeding 2,500 ng/mL. Therefore, false negative results may occur in users who have indicated or whose clinical status or history would indicate they are currently taking high doses of biotin (> 10 mg per day).
- False negative results are more likely after 4 days or more of symptoms.
- A negative test result may occur if the target antigen concentration in sample is below the detection limit of the

10. Rotate the swab 10 times.


- Do not lift the swab until you fully complete 10 rotations.

11. Remove the swab immediately from the sample port after rotating the swab.

- False negative results may occur if the order of testing procedures is not correctly followed.

12. Start a timer and read the result at 10 minutes. The test result should not be read after 15 minutes.


- Do not move or lift the test cassette. Results should not be read before 10 minutes or after 15 minutes.

PART D Interpretation of Results


- A control line next to "Cont" may appear in the result window in a few minutes but a test line may take as long as 10 minutes to appear. Incorrect results may occur if the test result is read before 10 minutes or after 15 minutes.

Positive Result


Influenza A Positive

One purple-colored line next to "Cont" and one red-colored line next to "A" indicate Influenza A positive result (Influenza A antigen detected)

Influenza B Positive

One purple-colored line next to "Cont" and one blue-colored line next to "B" indicate Influenza B positive result (Influenza B antigen detected)


COVID-19 Positive

One purple-colored line next to "Cont" and one blue-colored line next to "C19" indicate COVID-19 positive result (COVID-19 antigen detected)


- Look very closely! The color intensity in the test region will vary. Any faint colored line in the test region should be considered as positive.

Co-infection Positive Result

Any lines next to "A", "B", and/or "C19" with a control line next to "Cont" indicate Influenza A, B, and/or COVID-19 positive results. (Influenza A, B, and/or COVID-19 antigen detected)


Negative Result

One purple-colored line ONLY next to "Cont" indicates a negative result, with no test line next to "A", "B", and "C19". This means Influenza A, B, and COVID-19 antigen have not been detected.

- Negative results are presumptive and may need to be confirmed with a molecular assay.


Invalid

No purple-colored line next to "Cont" indicates an invalid result. A new test is needed to get a valid result. Repeat the test with a new sample and new test kit materials.

PART E Disposal

Dispose of all used test kit components and sample swab in household trash.

- test or if the sample is collected or handled improperly.
- False positive test results are more likely when the prevalence of SARS-CoV-2, influenza A and/or influenza B is low in the community.
- Viral Transport Media (VTM) should not be used with this test.
- This test does not distinguish between SARS-CoV and SARS-CoV-2.
- Persons with risk factors for severe disease from respiratory pathogens (e.g., chronic lung or heart disease, compromised immune system, diabetes, and other conditions) should contact a healthcare provider.
- This test is read visually. Because test lines can be very faint, users with conditions affecting their vision - such as far-sightedness, glaucoma, or color blindness - are encouraged to seek assistance to interpret results accurately (e.g., reading glasses, additional light source, or another person). This test has not been validated for use by those with color-impaired vision.
- This device is a qualitative test and cannot provide information on the amount of virus present in the specimen.
- This test detects both viable (live) and non-viable influenza A, influenza B, and SARS-CoV-2. Test performance depends on the amount of virus (antigens) in the sample and may or may not correlate with viral culture results performed on the sample.
- This test was not evaluated for use with the FluMist® vaccine. Individuals who recently received nasally administered influenza A or influenza B vaccine may have false positive test results after vaccination.

Frequently Asked Questions

- What are the known and potential risks and benefits of the test?**
Potential risks include:
 - Possible discomfort during sample collection.
 - Possible incorrect test result (see Warnings and Result Interpretation sections for more information).Potential benefits include:
 - The results, along with other information, can help you and your healthcare provider make informed recommendations about your care.
 - The results of this test may help limit the spread of COVID-19 and flu to the family of the tested individual and others in your community.
- What is the difference between an antigen and molecular test?**
There are different kinds of tests for the viruses that cause COVID-19 and the flu. Molecular tests detect genetic material from the virus(es) whereas antigen tests, such as the CareSuperb™ COVID-19/Flu A&B

Antigen Combo Home Test, detect proteins from the virus(es). Due to the lower sensitivity of antigen tests, there is a higher chance this test will give you a false negative result when you have COVID-19, Flu A and/or Flu B than a molecular test would.

- What if I have a positive test result?**
A positive result means that it is very likely you have COVID-19, Flu A and/or Flu B because proteins from the virus that causes COVID-19, Flu A and/or Flu B were detected in your sample. You should self-isolate from others and contact a healthcare provider for medical advice about your positive result.
 - What if I have a negative test result but still have symptoms?**
A negative test result indicates that antigens from the virus that causes COVID-19, Flu A and/or Flu B were not detected in your sample. If you have a negative result, it does not rule out SARS-CoV-2, Flu A and/or Flu B infection; you may still be infected and you may still infect others. It is important that you work with your healthcare provider to help you understand the next steps you should take.
 - What does an invalid test result mean?**
An invalid test result means the test did not work properly. If the test is invalid, a new swab should be used to collect a new nasal specimen and test again with a new test.
 - How accurate is this test?**
For more information on the performance of the test and how the performance may apply to you, please refer to the performance data in the Healthcare Provider Instructions for Use (IFU), available at www.accessbio.net/products/rapid-diagnostic-tests/caresuperb-covid-19-and-flu-a-and-b-antigen-combo-home-test.
 - Do I have to report my results? (optional)**
Report your test result(s) at [MakeMyTestCount.org](https://report.makemytestcount.org) - this voluntary and anonymous reporting helps public health teams understand COVID-19 spread in your area and across the country and informs public health decisions. <https://report.makemytestcount.org>.
 - What if I am uncertain how to interpret my results?**
If uncertain how to proceed, contact Technical Assistance at +1-888-898-1270 (Toll Free) (8 a.m. to 5 p.m. EST, Monday through Friday) or TShelp@accessbio.net.
- Important: Do not use this test as the only guide to manage your illness. Consult your healthcare provider if your symptoms persist or become more severe. Individuals should provide all results obtained with this product to their healthcare providers.

Scan the QR code to access the instructional video for a step-by-step guide on how to use this product.

Explanation of Symbols

- Single use
- Use by
- Manufacturer
- Date of manufacture
- Batch code
- Catalogue number
- Over-the-Counter
- Temperature limitation
- Contains sufficient for <n> tests
- Do not use if package is damaged
- Consult instructions for use
- In vitro diagnostic medical device

Manufacturer
Access Bio, Inc.
7 Fitzgerald Ave., Monroe Township, NJ 08831, USA
Tel: 732-873-4040 Fax: 732-873-4043 Web: www.accessbio.net
Technical Support
Tel: 888-898-1270 (Toll Free) E-mail: TShelp@accessbio.net

Quick Reference Instructions

CareSuperb™ COVID-19/Flu A&B Antigen Combo Home Test

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